



Mercury Thermometer

Weiyou® Mercury Thermometer is made of high quality medical glass and mercury. Weiyou® Mercury Thermometer is mainly used for Medical, Clinical offices, Pharmacy, Laboratory, Hospital, etc. As the manufacturer in China, Weiyou® has mature production lines and supply chains, and always supply best quality product, best favourable price and best service.

Product Description

Weiyou® Mercury Thermometer is made of high quality medical glass and mercury. Weiyou® Mercury Thermometer is mainly used for Medical, Clinical offices, Pharmacy, Laboratory, Hospital, etc. OEM customization is acceptable. Free samples, and small quantity can be supplied.

Specification

Item	Mercury Thermometer
Usage	Body or Animal Temperature Measurement
Material	Medical glass and mercury
MOQ	500pcs
Use For	Medical, Clinical offices, Pharmacy, Laboratory, Hospital, etc.
Advantage	High Sensitivity
Color	White, yellow, blue and other colors
Packing	Carton

Features

Material: Mercury liquid

Temperature range: 35-42 in C scale / 94-108 in F scale

Description: Clinical Thermometer for Armpit use





www.nbweiyou.com



Ningbo Weiyou Import & Export Co., Ltd.
Phone:+86-13567833512
E-mail:_dario@nbweiyou.com



www.nbweiyou.com



Ningbo Weiyou Import & Export Co., Ltd.
Phone: +86-13567833512
E-mail: dario@nbweiyou.com





www.nbweiyou.com



Ningbo Weiyou Import & Export Co., Ltd.
Phone: +86-13567833512
E-mail: dario@nbweiyou.com



